## Warranty

The manufacturer provides warranty in accordance with the legislation of the customer's own country of residence, with a minimum of 12 months, starting from the date on which the appliance is sold to the end user.

The warranty only covers defects in material or workmanship, which will not apply in cases of:

- Normal wear and tear
- Incorrect use, e.g. overloading of the appliance, use of non-approved accessories
- Use of force, damage caused by external influences
- Damage caused by non-observance of the user manual, e.g. connection to an unsuitable mains supply or non-compliance with the installation instructions
- Partially or completely dismantled appliances.

## Declaration of Conformity

This product is in accordance with the following European directives:

2014/30/EU (EMC)

2014/53/EU (RED)

2011/65/EU (RoHS)

## Disposal Information

If the icon below (crossed-out trashcan on wheels) is found on the product, the European directive 2002/96/EG applies. These products may not be disposed of in the household waste at the end of its useful life. Disposal can take place at appropriate collection points provided in your country.

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# SMART SCALE

## USER MANUAL

Read this manual thoroughly before use and save it for future reference

Fitleap  
Smart Scale

Find help and advice at  
[scale.fitleap.de](https://scale.fitleap.de)  
or scan the QR code on the left

## Technical Data

Model	Fitleap Smart Scale
Rated Voltage	4 x 1.5V AAA
Rated Power	0.6 W
Article No.	63866

### Note

Please contact the Fitleap customer support team if you have any questions or problems. You can reach us via [hello@fitleap.de](mailto:hello@fitleap.de). We're happy to help!

## Features

Thank you for purchasing the Fitleap Body Fat Scale. For your own safety, please read this manual carefully BEFORE using the device. Using the body fat scale has many advantages:

- Precise measurement of weight, body fat, muscle mass, water percentage, pulse and many other values.
- Management of up to 24 different user profiles via the associated app.
- Bluetooth synchronization with the app, so that weight histories and detailed information can be conveniently displayed and exported
- Large display with indication of the most important values

## Important safety precautions

### Warning

1. Do not use the device as a toy. Special attention is required when using near children.
2. Never immerse the device in water or other liquids. Never pour water over the device.
3. Avoid sources of heat and direct sunlight to prevent the device from heating up excessively.
4. Do not use the device for any purpose other than personal weighing.

## Requirements for using the scale

### Internet-enabled iOS or Android smartphone

For the following operations, you need a smartphone with iOS or Android operating system. A WLAN or 3G/4G/5G Internet connection and active Bluetooth are required to perform the following functions:

- Setting up the Fitleap scale
- Measurement of bio-measurement data with the Fitleap scale
- Retrieving / synchronizing measurement data from the Fitleap scale in general
- Access to certain functions of the Fitleap scale

### Smartphone operating system

The mobile app associated with Fitleap is the **Fitdays app**, which can be downloaded from the Apple App Store or the Google Play Store.

To use the Fitdays app on an iOS device, at least the iOS 8 operating system must be installed. The use of Bluetooth must be allowed.

To use the Fitdays app on an Android device, at least the Android 6.0 operating system must be installed. The use of Bluetooth and location services must be allowed.

We generally recommend keeping the operating system and the Fitdays app as up-to-date or fairly up-to-date as possible for optimal performance.

### Download Fitdays App

We recommend that you always install the latest version of the Fitdays app to take advantage of all the features. You can find the app in the Apple App Store or Google Play Store under "Fitdays" or by scanning one of the following QR codes.


Fitdays for iOS

Fitdays for Android


## Product diagram

1 Display

2 Electrodes


## Battery compartment

The battery compartment is located on the back of the scale and can be opened by sliding it open, applying a little pressure on the battery compartment lid if necessary.

To operate the scale, 4 AAA batteries are required, which are inserted in alternating direction according to the marking in the battery compartment.

## Display diagram

- 1 Battery indicator
- 2 Bluetooth connection status
- 3 Weight unit
- 4 Weight
- 5 Weight distribution left/right
- 6 Body fat percentage

It is possible that the letters below can be seen on the display. They have the following meaning:

Display	Meaning
(Display dark)	The scale has no power. Replace the batteries.
Lo	The battery power is weak or decreasing. Replace the batteries.
EEr	The scale is overloaded. Reduce the weight.
[	The balance is calibrated. Wait a little while until the calibration is completed.

- 7 Water percentage
- 8 Pulse
- 9 Body Mass Index (BMI)
- 10 Muscle percentage
- 11 Bone percentage

## Before first use

Before using the scale, the scale must be positioned correctly in order to obtain accurate measurement results.

The next step is to connect the scale to the Fitdays app to activate and use all the functions of the scale.


### NOTE

All functions of the scale, such as body fat measurement, are only active and usable after pairing the scale and app and creating a user account. An anonymous guest account can also be used for this purpose.

## Preparations for precise weighing

### The correct floor

In order to provide correct weighing results, it is important that the scale stands securely on a level surface. Smooth, solid floors such as tiles or stone and wooden floors in general are suitable for this purpose. The scale should not stand on a carpet, be tilted or wobble.

### Standing correctly

When weighing, you should stand on the scale barefoot and with dry feet, otherwise the body fat measurements, among other things, may deviate from the real values. Make sure that your feet are parallel to each other and centered in the measuring ranges of the scale.

## Set up the app and connect to the scale

First install the Fitdays app on your smartphone (see p. 16) and then follow the steps below to set up the app correctly and connect it to the scale.



To do this, first open the Fitdays app on your smartphone. When you first open the app, you will be asked to allow the use of Bluetooth. This is necessary for the correct connection with the scale.

## Connect scale and app


You can now optionally use the app either as a guest, or register a personal account or log in to an existing account. When registering an account, data such as gender, height and age are entered, which have an influence on the measured or calculated values. When using as a guest, your data is not saved, so if you change smartphones or reinstall the app, it will be lost unless it has been exported beforehand.

After entering your data, the scale can now be added as a device in the app. To do this, proceed as follows:

1. Tap "User account" on the bottom right of the main screen.
2. You are now in the profile overview. Tap on "Device" at the top left to switch to the overview of connected devices.

- On the following screen, tap the plus icon on the top right to start pairing your scale with the app.
- Tap the scale with your foot or stand lightly on the scale so that the display lights up. The scale is ready to be paired. Now select the Pairing with Bluetooth option by tapping the Bluetooth icon that appears on the screen in the Fitdays app. The scale will automatically connect to the app and will be displayed in the device overview after successful pairing.

**Note**

If the connection between the app and the scale via Bluetooth does not work, the app will display a message after a short time. If this does not solve the problem, please contact the FitLeap Customer Support Team at [hello@fitleap.de](mailto:hello@fitleap.de). We will be happy to help you!

## Weighing with the scale

After you have performed the steps described above, the scale is ready for weighing. Stand barefoot on the scale as described on page 19. The scale starts measuring.

### Even weight distribution

The weight distribution indicator shows you if you are standing balanced or if you are shifting your weight too much to one side. If this is the case, try to balance it out.

The weight distribution indicator is visible both on the scale and on the app.

The weight distribution indicator is visible both on the scale and on the app.

### Body fat, pulse and other values


Following the weight measurement, the scale records and calculates other values such as pulse, body fat, BMI, muscle mass and bone mass. To do this, stay on the scale until the measurement is finished and all values are recorded in the app.

The results of the measurement are displayed in the lower half of the screen, where a detailed view can be accessed by tapping the "More" button at the bottom.

### Detailed history

In the app, from the main menu, tap "Graph" in the bar at the bottom of the screen to view detailed progressions of various readings over time.

You can also activate the week, month or year view by tapping the corresponding selection above the graph.



## Multiple user profiles

Through the Fitdays app, up to 24 different user profiles can be connected to the Fitleap body fat scale at the same time.

### Adding users and switching between users


To add a new user profile, select "User Account" in the bottom right of the main screen, then "Add User" in the top center and enter the appropriate data. By "swiping left" on a profile in the list, it can be deleted again.

Switching between different user profiles is done via the main menu by tapping on the currently active profile name at the top left of the screen and selecting the desired user in the menu that pops up.

### Assigning weighing results


If a weighing has taken place which cannot be clearly assigned to a user profile (usually due to the large weight deviation), then it is displayed in the app and can be assigned manually.

## The baby mode

The Fitleap body fat scale, in combination with the Fitdays app, allows you to weigh babies when they are being carried by adults. To weigh a baby, select the "+" at the top right of the main screen and then select "Baby Mode".

The next screen is where the weighing takes place: First, the adult person weighs himself alone and then with baby in his arms. The scale will automatically calculate the corresponding weight.

## Recalibrate the scale

If you have moved the scale to a new surface or position, or turned it over (e.g., to change the batteries), it is almost always recommended that you recalibrate it before weighing. You calibrate the scale as follows:

1. Place one foot on the scale and apply a little pressure so that the display lights up. Then remove the foot from the scale again.
2. The balance recalibrates itself. When the display shows 0.00, the calibration is complete and you can use the scale as usual.

## Import data to Apple Health, Google Fit and fitbit

You have the option to import your measurement results from the Fitdays app into Apple Health, Google Fit and fitbit.

For this purpose, we provide you with step-by-step instructions, which you can view online for convenience by visiting the corresponding webpage below or by scanning the corresponding QR code.

If you have any problems connecting, please feel free to contact our Customer Support Team at [hello@fitleap.de](mailto:hello@fitleap.de).

You can find the instructions on

Apple Health

<https://apple.fitleap.de>

Google Fit

<https://google.fitleap.de>

Fitbit

<https://fitbit.fitleap.de>

## Problem Solving

Problem	Possible Cause	Solution
The scale cannot be connected via Bluetooth.	The operating system of your smartphone is not up to date.	Update the operating system to a correspondingly current version
	The device is damaged.	Contact Fitleap Customer Support
Fat is not measured during weighing.	You do not stand on the scale barefoot and dry feet	Make sure you are barefoot and both the scale and your feet are dry.
	Weighing is not finished yet	Wait until the weighing process is displayed as completed.
	Bluetooth is not enabled	Activate Bluetooth on your smartphone and open the Fitdays app.
	The app is paired with another scale/device.	Unpair the other scale and pair the app with the Fitleap scale again.

## Problembehebung

Problem	Mögliche Ursache	Lösung
The scale shows an inconsistent weight.	The scale sits on a carpet or soft surface.	Check that the scale is standing on a hard, level surface.
	Dirt or material is stuck to the sensor feet that will interfere with accurate measurement.	Check each sensor foot on the back to make sure nothing is stuck to the bottom. Then recalibrate the scale.
The app shows the wrong fat to muscle ratio.	The setting of the body parameters is incorrect.	Check that you have entered the correct gender, size and age.
The scale will not turn on.	The batteries are not positioned correctly or are not seated correctly.	Remove the batteries from the scale and reinsert them. Be sure to insert the batteries in the correct direction.
I would like to recalibrate the scale.		Step on the scale to turn it on. Automatically drop the scale to 00 kg to calibrate it.
The measured values for water and body fat percentage vary too much.	Weighing at different times of the day can affect body fat percentage readings.	Weigh yourself at the same time every day. The body fat percentage is always higher in the morning. The water percentage is always higher in the evening.
The scale shows only weight, balance and BMI.	You may not weigh yourself barefoot or with dry feet / on a dry scale.	Weigh yourself barefoot and with dry feet / on dry scales.
	You may not have your own user account in the Fitdays app or the app may not be paired with the scale	Make sure you have your own account in the app and have weighed yourself at least once with the app connected before weighing yourself offline.

### Note

Please contact the Fitleap Customer Support Team at [hello@fitleap.de](mailto:hello@fitleap.de) with any questions or concerns. We are happy to help!