Warranty

The manufacturer provides warranty in accordance with the legislation of the customer's own country of residence, with a minimum of 12 months, starting from the date on which the appliance is sold to the end user.

The warranty only covers defects in material or workmanship, which will not apply in cases of:

- Normal wear and tear
- Incorrect use, e.g. overloading of the appliance, use of non-approved accessories
- Use of force, damage caused by external influences
- Damage caused by non-observance of the user manual, e.g. connection to an unsuitable mains supply or non-compliance with the installation instructions
- Partially or completely dismantled appliances.

Declaration of Conformity

This product is in accordance with the following European directives:

2014/30/EU (EMV)

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Disposal Information

If the icon below (crossed-out trashcan on wheels) is found on the product, the European directive 2002/96/EG applies. These products may not be disposed of in the household waste at the end of its useful life. Disposal can take place at appropriate collection points provided in your country.

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
MASSAGE GUN

USER MANUAL

Read this manual thoroughly before using and save it for future reference

Fitleap
Massage Gun Pro

Find help and advice at
massage.fitleap.de
or scan the QR code on the left

Technical details

Model	Fitleap Massage Gun Pro
Rated Voltage	220-240 V / 50-60 Hz
Rated Power	25 W
Article No.	72476

Note

Please contact the Fitleap customer support team if you have any questions or problems. You can reach us via hello@fitleap.de. We're happy to help!

Features

Thank you for purchasing the Fitleap Massage Gun. For your own safety, please read this manual carefully BEFORE using or charging the device. There are many benefits to using this massage gun:

- Tense or fatigued muscles can be loosened during exercise
- In case of muscle soreness after high exertion, the massage gun can be used to massage lactate out of the muscle tissue
- The massage gun can also be used outside of sports to relieve severe tension
- Different attachments and massage intensities for specific muscle groups help to find exactly the right combination for a particular muscle.

Important safety precautions

Warning

Before using the massage gun on your body, please get your doctor's approval first if any of the following factors apply to you:

- Pregnancy, diabetes with complications (e.g., neuropathy or retinitis), wearing a pacemaker
- Peripheral neuropathy, recent surgery, epilepsy, migraines
- Herniated disc, spondylolisthesis, spondylolysis, or spondylosis,
- Recent joint replacement or IUDs, metal pins or plates, or other concerns about your physical health and constitution.

In general, the fitness gun is not suitable for use by children or anyone under the age of 13.

Charging the battery

1. Make sure that the device is switched off. The battery indicator should not be lit. If this is not the case, switch off the device by long pressing the POWER button.
2. Plug the charger into a suitable socket and then connect the charging cable to the charging socket provided under the handle of the massage gun.
3. When the cable is correctly connected, the charging process begins, indicated by the charging indicator in the handle, which constantly changes color. A full charge of the battery usually takes between 2 and 3 hours.
4. When the battery is fully charged, the charging indicator in the handle does not change color any further, but lights up green continuously. Now remove the charging cable from the massage gun. The charging process is complete and the device can be used.

In general, the battery can be charged in any state of charge, not only when completely discharged.

Battery charge states


Battery indicator color	Meaning	Solution
The charging indicator flashes red	The battery of the device is almost completely empty. The device will soon switch off automatically.	Stop operation, turn off the device and fully charge it via the charging cable.
The charging indicator lights up red	The battery of the device is running low.	Fully charge the device promptly via the charging cable.
The charging indicator lights up in changing colors	The battery of the device is being charged. The color change indicates the active charging process.	Allow the device to continue charging. When the indicator turns solid green, charging is complete.
The charging indicator lights up permanently green	The battery is fully charged.	Disconnect the charging cable. The battery is now fully charged and the device is ready for use.

IMPORTANT

For optimal battery life, never fully discharge the device until the battery is completely empty. Instead, it is best to charge it immediately when the charge indicator turns red and then fully charge the battery.

The different massage modes


By pressing the POWER button several times, you can switch between 6 different massage modes. The active mode is indicated by a brightly lit symbol.

Display Icon	Title & Speed	Meaning
	Smart 1800 - 2500	In Smart mode, the speed of the massage varies in order to stimulate the muscle in different ways.
	Training 1650	The low speed training mode is designed for use during workouts.
	Relax 1900	Relax mode is designed for the recovery phase after training.
	Relief 2100	Relief mode is designed to help with sore muscles by massaging lactates out of the muscles.
	SPA 2300	SPA mode is suitable for general massages, regardless of training or sports.
	Fitness 3200	Fitness mode has the highest speed and is especially designed for athletes whose muscles are used to high load.

Using the massage gun

1. Select the desired massage attachment and insert it with a little pressure into the recess on the front of the massage gun.
2. Turn on the massage gun by long pressing the POWER button. The massage gun is ready for use when the battery indicator in the handle is not illuminated red (an indicator illuminated blue indicates that the battery is sufficiently charged).
3. Decide on a massage mode and press the POWER button on the display until the desired mode is active. The respective active mode is indicated by the white lit symbol.

If you are not sure which mode is best for the muscle group you want to massage, select Smart mode. In this mode, the massage speed varies in order to stimulate the muscle in different ways.

WARNING

Do not operate the massage gun if it is currently being charged via the power supply unit. Ideally, wait until the battery is fully charged before disconnecting the device from the charger and using it.

The intensity and duration of the application

The intensity and duration of use of the massage gun depends primarily on the phase of training in which it is used. In general, one can orientate oneself on the following areas here

Muscle Activation *(Duration: max. 30 seconds per muscle group)*

Before a workout or activity, muscles and muscle groups can be worked for up to 30 seconds using the massage gun. This activates muscles and nerves.


Muscle Re-Activation *(Duration: max. 15 seconds per muscle group)*

If your muscles are already feeling a bit tired in a workout, and the muscles are contracting or cramping, you can try to loosen them back up a bit with a short reactivation.

Recovery and pain relief *(Duration: max. 2 minutes per muscle group)*

After the workout or when the soreness is already noticeable, the massage gun can be used for up to two minutes per muscle group to loosen the muscle tissue. This can be repeated several times a day, with at least 1-2 hours between applications.


The massage attachments

Medium

About the Medium attachment


Made of slightly yielding high-density foam, this attachment is designed for general use, but can also be used around bones and other delicate areas.

Bullet

About the Bullet attachment

The small diameter of the bullet makes it perfect for targeted muscle treatment, especially in deeper muscle regions. Ideal for localizing and eliminating trigger points.


Fork

About the Fork attachment

Optimally suited both for the treatment of muscles around tendons and ligaments, and for therapy around the spine, shoulder and calf.


The massage attachments

Wide

About the Wide attachment


The Wide Head is designed to stretch out the muscles. By working in one direction, this attachment can "push" lymph fluid and lactate out of the muscles.

Hammer

About the Hammer attachment

The hammer attachment can be used to work larger, less sensitive muscle groups such as thighs, which can withstand a little more pressure. Likewise, it is well suited for warm-up and activation

Soft

About the Soft attachment

The Soft attachment is particularly suitable for more sensitive areas due to its soft shape, so it should be used near tendons or smaller muscle groups, for example.

Cleaning and storage

Never clean the device while it is charging and make sure that the device is switched off. This will prevent the risk of electric shock.

Do not use water or harsh or corrosive cleaners or solvents for cleaning. Instead, take a slightly damp cloth and gently wipe the individual components of the massage gun (with the exception of the charger and the charging socket).

Transport

If possible, the massage gun should not be transported over long distances fully loaded.

Problem solving

Problem	Possible cause	Solution
The massage gun does not turn on	The battery is empty	Laden Sie den Akku über das Netzteil vollständig auf und versuchen Sie es erneut
	The device is damaged	Contact the Fitleap Customer Service team
The massage gun works very slowly / very sluggishly	The battery is almost empty	Fully charge the battery via the power supply unit and try again
	The battery capacity has decreased over the years	Contact the Fitleap Customer Service team
	The device is damaged	Contact the Fitleap Customer Service team
The massage gun can no longer be charged	The charging cable is not fully inserted into the charging socket	Pull out the charging cable and insert it completely into the socket
	The charger, the charging cable or the battery are defective	Contact the Fitleap Customer Service team
The display has scratches	The protective cover foil has not been removed yet	Remove the protective cover foil from the display